

700 Series MICROLITER Syringes

5 μ L - 500 μ L

- Cemented needles (N), removable needles (RN), or luer tip (LT)
- Plungers and syringe barrels are not interchangeable or replaceable
- Tight tolerances between the plunger and the barrel



- For use with liquids
- All syringes are autoclavable when disassembled except for the cemented needle (N) syringes











700 Series MICROLITER Six Pack Syringes

10 μ L Six Packs

- Savings! Syringes cost less in six packs
- Convenience! Packaging container lets you keep syringes close at hand
- Identification! Six numbered Ident-A-Dots provide quick syringe identification



Ordering Information

		Volume Model	5 µL	10 µL	25 µL	50 µL	100 µL	250 µL	500 µL
			75	701	702	705	710	725	750
			Gauge**	26s	26s	22s	22s	22s	22s
N (Cemented Needle)			87900	80300	80400	80500	80600	80700	80800
N (Cemented Needle)			87919	80383	80465	80565	80665	80765	80865
SN (Cemented Needle) (Specify needle gauge, length, and point style)			87908	80308	80408	80508	80608	80708	80808
RN (Removable Needle)			87930	80330	80430	80530	80630	80730	80830
RN (Removable Needle) Syringe (Needle not included)			7634-01	7635-01	7636-01	7637-01	7638-01	7639-01	7640-01
LT (Luer Tip) (Needle not included)				80301	80401	80501	80601	80701	80801
PB600 Repeating Dispenser (Syringe not included)					83700	83700	83700	83700	83700
Reproducibility (Chaney) Adapter (Syringe not included)			14700	14700	14725	14725	14725	14725	14725
Syringe Guide (Syringe not included)			14806	14806	14906	14906	14906	14906	14906
Digital Base Unit (Syringe not included)			Add the prefix "DS" to the beginning of the syringe's part number. (Base unit not sold separately)						



Syringe Calibration Service


Add the prefix "CAL" to the beginning of the syringe's part number.

Typical Operating Parameters



- Maximum temperature for Cemented Needle 50° C (122° F) and all others 115°C (239° F).
- Maximum pressure rating is 2000 psig for 5 µL to 100 µL, all others are 1000 psig.

Replacement Parts

Description	Volume	5 µL	10 µL	25 µL	50 µL	100 µL	250 µL	500 µL
	Model	75	701	702	705	710	725	750
1. RN Nut		30902	30902	30902	30902	30902	30902	30902
2. RN Needle, see page 31								



Ordering Information

		Volume	10 µL
		Model	701
		Gauge**	26s
N (Cemented Needle)			80366
RN (Removable Needle)			80336