

Syringe Accessories

MTB Sampling Valve

- Lever-actuated on/off valve for storing and transporting samples in any GASTIGHT TLL syringe
- Kel-F hub needles attach to the valve
- Pressure tested to 100 psig; autoclavable
- Bore I.D., .039"



Septum Adapter

- Store and sample syringe contents through the Septum Adapter
- Attaches to any TLL syringe or luer lock fitting
- Replaceable septum



Removable (RN) Adapters

- For 250 μ L and larger removable needle (RN) syringes
- Port I.D. of .039" allows sample spiking with a 701N syringe
- Provides connection to threaded fittings, tubing assemblies, tees and crosses



Ordering Information

Part Number	Description
86580	MTB Sampling Valve

Ordering Information

Part Number	Description
31335	Septum Adapter (TLL), includes 12 septa
75810	Replacement Septa, 7 mm dia., 12/pk

Ordering Information

Part Number	Description
35081	Female Luer/Removable Needle Hub (L)
35080	Male Luer/Removable Needle Hub (L)
35083	Male Luer Lock/Removable Needle Hub (L)