

Syringe Graduations

Syringe Graduations

100 line screens

7000 Series MODIFIED MICROLITER Syringes

	.5 μ L	1 μ L	2 μ L	5 μ L
	7000.5	7001, 7101	7002, 7102	7105
Scale length, cm	6	6	6	6
Total scale divisions	100	100	100	100
Major graduations, μ L	0.05	0.1	0.2	0.5
Minor graduations, μ L	0.01	0.02	0.04	0.1
Sub-minor graduations, μ L	0.005	0.01	0.02	0.05

600, 700, 800 & 900 Series MICROLITER Syringes

	2.5 μ L	5 μ L	5 μ L	10 μ L	25 μ L	50 μ L	100 μ L	250 μ L	500 μ L
	62	65	75, 85, 95	701, 801, 901	702, 802	705, 805	710, 810	725, 825	750
Scale length, cm	2.7	2.7	5.41	5.41	6	6	6	6	6
Total scale divisions	50	50	100	100	100	100	100	100	100
Major graduations, μ L	0.5	1	0.5	1	2.5	5	10	25	50
Minor graduations, μ L	0.1	0.2	0.1	0.2	0.5	1	2	5	10
Sub-minor graduations, μ L	0.05	0.1	0.05	0.1	0.25	0.5	1	2.5	5

1700 & 1800 Series GASTIGHT Syringes

	10 μ L	25 μ L	50 μ L	100 μ L	250 μ L	500 μ L
	1701, 1801	1702, 1802	1705, 1805	1710, 1810	1725, 1825	1750
Scale length, cm	6	6	6	6	6	6
Total scale divisions	100	100	100	100	100	100
Major graduations, μ L	1	2.5	5	10	25	50
Minor graduations, μ L	0.2	0.5	1	2	5	10
Sub-minor graduations, μ L	0.1	0.25	0.5	1	2.5	5

1000 Series GASTIGHT Syringes

	1 mL	1.25 mL	2.5 mL	5 mL	10 mL	25 mL	50 mL	100 mL
	1001	1001.25	1002	1005	1010	1025	1050	1100
Scale length, cm	6	6	6	6	6	6	6	12
Total scale divisions	100	100	100	100	100	100	100	200
Major graduations, mL	0.1	0.125	0.25	0.5	1	2.5	5	5
Minor graduations, mL	0.02	0.025	0.05	0.1	0.2	0.5	1	1
Sub-minor graduations, mL	0.01	0.0125	0.025	0.05	0.1	0.25	0.5	0.5

Super Syringes

	500 mL	1000 mL	1500 mL	2000 mL
	S0500	S1000	S1500	S2000
Scale length, in.	6.53	6.57	9.85	13.13
Total scale divisions	25	50	75	100
Major graduations, mL	100	100	100	100
Minor graduations, mL	20	20	20	20